

Clinical Experience with the Xchange™ Nasal Dock in Deep Sedation for GI Procedures

Stephen Garber, MD

Introduction

Airway obstruction leading to oxygen desaturation and hypercapnia remains one of the most common and serious complications encountered during moderate to deep sedation—particularly in patients with preexisting airway pathology such as obesity, obstructive sleep apnea, or anatomic abnormalities. Maintaining airway patency and adequate gas exchange during such procedures often requires repeated interventions, additional personnel, and occasionally conversion to general anesthesia.

The Xchange™ Nasal Dock (Pneuma Therapeutics) was developed to address this problem directly. It is a single-use nasal CPAP device engineered to provide rapid, effective continuous positive airway pressure (CPAP) throughout both the intraprocedural and post-procedural periods. The system attaches easily to standard wall oxygen, incorporates a simple adjustable APL valve for pressure titration, and provides visual feedback of ventilation via a reservoir bag.

Clinical Experience

In an initial clinical evaluation involving 10 patients undergoing deep sedation for gastrointestinal procedures, all with varying degrees of known airway pathology, the Xchange™ Nasal Dock demonstrated consistent effectiveness and ease of use.

Key Findings

- 1. Ease of Application and Rapid Effectiveness:** Application was intuitive and required minimal setup time. Effective positive pressure was achieved within 20 seconds of placement.
- 2. Patient Comfort and Acceptance:** Patients tolerated the device well, reporting minimal discomfort compared to traditional nasal CPAP masks or oral airways.
- 3. Efficient Oxygen Source Utilization:** The device operated effectively using standard wall oxygen, eliminating the need for external flow regulators or specialized equipment.
- 4. Simple Pressure Adjustment:** The integrated APL valve allowed straightforward titration of airway pressure, facilitating rapid response to clinical needs.
- 5. Visual Verification of Ventilation:** The attached anesthesia reservoir bag provided excellent visual confirmation of ventilation and respiratory effort throughout the procedure.

6. Stable Oxygenation and Ventilation: Across all cases, oxygen saturation and end-tidal CO₂ levels were maintained within optimal ranges, even with patients positioned laterally—conditions that often challenge airway maintenance.

Summary and Conclusion

This early clinical experience suggests that the Xchange™ Nasal Dock is an innovative, practical, and effective solution for maintaining airway patency and adequate ventilation in patients receiving moderate to deep sedation. The device's rapid setup, simplicity, and visual feedback features enhance both provider confidence and patient safety.

While this report reflects a small initial sample, these findings indicate that the Xchange™ Nasal Dock will likely occupy an important role in the procedural sedation landscape, particularly for patients with known airway pathology where airway management remains a critical concern.